


# AgarSqueezer – a cell confiner


Confine your cells under agarose nanopillars to study cell behaviour



**AgarSqueezer** is a device designed to study cell response to short and long-term mechanical confinement within a physiological rigidity range.



Cell without confinement



Cell compressed within a thin space

Compression height:

- 2.5  $\mu\text{m}$
- 5  $\mu\text{m}$
- 30  $\mu\text{m}$
- 100  $\mu\text{m}$

## Key features

### 1. Instant & homogeneous cell confinement.

### 2. Physiological rigidity.

The mechanical properties of agarose can reproduce stiffness of the in vivo microenvironment (1-150 kPa).

### 3. Long-term confinement.

The porous nature of agarose facilitates nutrient and oxygen diffusion, allowing for long-term cell culture and monitoring in confined conditions (several days).


### 4. Fully compatible with in situ and ex situ analyses.

The system is compatible with real-time dynamic imaging and all immunostaining steps can be performed in situ. Alternatively, cells can be easily collected for standard molecular biology or functional assays.

### 5. Highly flexible.

Tunable pillar height, matrix stiffness & composition and possibility to coat with ECM proteins.

## Results from users



### Compression of immature TF1-GFP hematopoietic cells


Quantification of cell morphology under confinement. (A–C): Morphology of immature TF1-GFP hematopoietic cells for control (A) and for 30  $\mu\text{m}$  and 5  $\mu\text{m}$  (B and C, respectively). Scale bar = 20  $\mu\text{m}$ .

From A. Prunet et al. *Lab on Chip*, 2020.

### Arabidopsis root cells confined in Agarsqueezer

*Arabidopsis thaliana* Col-0 root cells stained with Calcofluor (cell wall) and imaged with a confocal microscope either in a traditional liquid culture (left), or after 24h of confinement under the 30 $\mu\text{m}$  (middle) or 5 $\mu\text{m}$  (right) pillars.

Image credits: Léa Bogdziewicz – UPSC – SLU Sveriges lantbruksuniversitet, Sweden



Original publication:

A. Prunet et al., A new agarose-based microsystem to investigate cell response to prolonged confinement, *Lab on a Chip*. 20:4016–4030 (2020)